

Monofocal  
Hydrophobic  
Preloaded



### Technical Specifications

Commercial name	MICROPURE 123		
Material	PhysIOL G-free® (GFY) (hydrophobic acrylic glistening-free) <sup>1</sup>		
Overall diameter	OD to 24.5D: 11.00 mm 25D to 30D: 10.75 mm		
Optic diameter	OD to 24.5D: 6.00 mm 25D to 30D: 5.75 mm		
Optic	Aspheric aberration-correcting (-0.11 $\mu$ SA)		
Filtration	UV & blue light		
Refractive index	1.52		
Abbe number	42		
Angulation	2°		
Injection system	PhysIOL 1.2.3		
Incision size	≥ 2.2 mm		
Spherical power	OD to 9D (1D steps) & 10D to 30D (0.5D steps) Cartridge with PRS® technology <sup>2</sup>		
Square edge	360°		
Nominal manufacturer A constant	119.40		
Suggested A constant <sup>3</sup>		<b>Interferometry</b>	<b>Ultrasound</b>
	<b>Hoffer Q: pACD</b>	5.85	5.59
	<b>Holladay 1: Sf</b>	2.06	1.80
	<b>Barrett: LF</b>	2.09	-
	<b>SRK/T: A</b>	119.40	119.05
	<b>Haigis<sup>4</sup>: a0; a1; a2</b>	1.70; 0.4; 0.1	1.214; 0.4; 0.1
	MICROPURE (non-preloaded)		
Overall diameter	-10D to 24.5D: 11.00 mm & 25D to 35D: 10.75 mm		
Optic diameter	-10D to 24.5D: 6.00 mm & 25D to 35D: 5.75 mm		
Injection system	Medical Accuject 1.8 up to 24.5D & Accuject 2.0/2.1/2.2 up to 35D		
Incision size	≥ 1.8 mm		
Spherical power	-10D to 9D (1D steps) & 10D to 30D (0.5D steps) & 31D to 35D (1D steps)		

<sup>1</sup> The PhysIOL G-free® (GFY) is patented since 2010. Chassain C, *J Fr Ophthalmol* 2018, 41(6):513-520.

<sup>2</sup> The PRS® technology is patent pending.

<sup>3</sup> Estimates only: surgeons are recommended to use their own values based upon their personal experience. Refer to our website for updates.

<sup>4</sup> Not optimized.

## Product Information


<b>Manufacturer</b>	PhysIOL s.a. - Liège Science Park Allée des Noisetiers 4 B-4031 Belgium +32 4 361 05 49 physiol@bvimedical.com
<b>Certificate information</b>	CE: Certificate N° CE658516 ISO 13485:2016: Certificate n° MD658518 MDSAP: Certificate N° MDSAP 691544 ISO 9001:2015: Certificate N° FM 658519
<b>Shelf life</b>	Three (3) years from manufacturing date for MICROPURE 123 Five (5) years from manufacturing date for MICROPURE
<b>Intended Use</b>	Intended use (for all IOLs): The posterior chamber intraocular lens which is intended to be placed into the capsular bag for the replacement of the human lens to achieve the visual correction of aphakia in adult patients in whom the cataractous lens has been removed by extracapsular cataract extraction.
<b>Indication for use</b>	The lens should be used as intended in adult patients surgically treated for cataract, with possibly associated presbyopia, who desire improved uncorrected far vision, with reduced spectacle dependence.
<b>Product Composition</b>	No products of animal or human origin are present in the implant. The implant is made of the GFY material proprietary to PhysIOL. It is composed of an acrylate copolymer Ethylene Glycol Phenyl Ether Acrylate (2-Phenoxyethyl Acrylate) (EGPEA) and 2 Hydroxyethyl Methacrylate (HEMA) including a UV light filter and a blue light filter
<b>For sterile product</b>	All IOLs from PhysIOL are steam sterilized
<b>Packaging Material</b>	Holder (Polypropylene) Container (Polypropylene) Storage liquid (0.9% NaCl solution) Aluminium lid (Aluminium Gold) Container label (paper) Blister PP (Polypropylene) Tyvek lid
<b>Product Class</b>	MDD Class IIb Sterile, According to European Medical Device Directive 93/42/EEC Not available in the United States

# MICROPURE 123 with Single Use Injector 1.2.3 Premium

## For 2.2-2.4 mm Incisions with PRS® (Pressure Release System) Technology

The MICROPURE 123 lens is delivered preloaded in a cartridge, which is simply clipped to the Single-Use Injector 1.2.3. Premium.

The Single-Use Injector 1.2.3. Premium requires no lens handling which ensures perfect control of asepsis and makes lens injection comfortable and reproducible. Additionally the unique PRS® technology offers an extremely smooth injection in combination with a significant decrease of pressure on the incision.



## Injection Guidelines

1. Connect the injector vertically onto the preloaded cartridge until you hear the “clip” indicating that both elements have been firmly and adequately locked. If you do not hear the “clip”, there is a possibility that the connection of both elements could not be secured. In the event you do not hear the “clip”, first remove the assembled device from the container. Secondly, vertically place the assembled device into the container and proceed once again to the “clipping”.
2. Push the plunger completely down towards the safety catch and...
3. ... keep the plunger in this position for 3 seconds. This ensures the lens is securely loaded in the cartridge. Then, gently release the plunger.
4. Remove the safety catch by a twist motion.
5. Remove the assembled system from the container and check that the cartridge is properly locked onto the injector. The non-return safety clip of the cartridge should be located just behind both lateral marks of the injector body, as illustrated in the above picture.
6. Rinse the IOL with Balanced Salt Solution (BSS) by introducing the cannula of the Balanced Salt Solution (BSS) syringe into the small hole on the body of the injector, and then inject a generous amount of ophthalmic viscoelastic device (OVD)<sup>1</sup> into the same hole.
7. Push the plunger for injection. When the first two haptics are out of the cartridge, release the plunger a few millimeters to free both posterior haptics, then push again until the implantation.


<sup>1</sup> Take the ophthalmic viscoelastic device (OVD) solution out of the refrigerator at least one hour before use.

## MICROPURE (non-preloaded) Injection System

The Medical Accuject 1.8 / 2.0 / 2.1 / 2.2 injection systems are recommended for implanting the MICROPURE (non-preloaded) lenses.

These fully single-use systems represent total reliability for safe and effective lens injections.

Their compact design with integrated cartridge enables a simple, predictable loading and positioning of the lens.



1. Apply ophthalmic viscoelastic device (OVD) into the tip and the loading chamber of the injector cartridge.
2. Remove the lens from the lens holder. Position the lens into the cartridge in such way that the two haptics with the notches are pointing at 1 and 7 o'clock.
3. Exert slight pressure onto the lens optic and make sure that all haptics are inside before further closing the cartridge. Close the cartridge and check the position of the lens.
4. Once the "click-lock" mechanism engages, the lens is securely loaded and ready for injection.
5. Press the injector plunger forward and push the lens into the conical tip of the cartridge.
6. Pull the plunger back a few millimeters and then inject the lens in one continuous motion. For gentle implantation, it is not necessary to fully push the plunger to the bottom of the cartridge.